


THE CLEARLY SMARTER

WAY TO STAY AHEAD OF HIGHS & LOWS

GUARDIAN™ CONNECT

The Continuous Glucose
Monitoring System for people on
multiple daily insulin injections



Medtronic
Further, Together

For a list of compatible operating systems and devices please visit: guardianconnect.medtronic-diabetes.co.uk/app-compatibility

Managing diabetes can feel like a guessing game. Will I go low on my drive to work? Am I too high after lunch? Are my levels falling too fast before I go to sleep? These questions you may ask yourself every day.

The Guardian™ Connect Continuous Glucose Monitoring (CGM) system is right by your side - helping you see what's really going on with your glucose levels.


CLEARLY SMARTER MONITORING HELPS YOU OUTSMART HIGHS & LOWS

Knowing what's ahead can change everything.

Now, there's a smarter continuous glucose monitor that can predict your highs and lows and alert you **before** they happen.

When you have the power to *predict* your glucose levels, you have the power to *take action*. So you can confidently focus on your life, not on your levels.


Guardian™ Connect also links seamlessly to CareLink™ Connect allowing you to share your data with up to 5 people who can also receive personalised alerts for added reassurance.



CLEARLY SMARTER GLUCOSE READINGS, RIGHT ON YOUR SMARTPHONE.

Now, you can check your glucose as easily as you check your text messages.


Guardian™ Connect continuously measures your glucose levels, and wirelessly sends readings to your smartphone using Bluetooth technology every 5 minutes (that's 288 times per day!).



CLEARLY SMARTER ALERTS LET YOU TAKE ACTION BEFORE A HIGH OR LOW.

Guardian™ Connect uses smart technology to analyse your glucose patterns and predict when you're going too high or low.


It's the **only** CGM that can alert you about predicted highs and lows anytime from 10 to 60 minutes in advance. So you can stop them, **before** they ever happen.



FUNDING OPTIONS

Medtronic are working with the NHS and HSE to improve patient access to CGM. Speak to your healthcare professional to see if you meet your local health authorities funding criteria.

For patients looking to self-fund Guardian™ Connect please visit our website for further information on our range of subscription packages.



CLEARLY SMARTER CGM TECHNOLOGY THAT FITS SEAMLESSLY INTO YOUR LIFE

The small and discreet sensor and transmitter can go with you anywhere – even swimming or in the shower, whilst the easy-to-use app lets you quickly see your glucose levels and receive customised alerts, keeping you in control.

GUARDIAN™ CONNECT FEATURES

- Insulin dose & blood glucose event markers
- Meal / snack tracking
- Activity and exercise tracking
- Add notes to your trends
- Personalised high and low alerts
- Set up to 5 people to receive personalised SMS alerts

CLEARLY CONNECTED TO YOUR CARE PARTNERS

You can stay connected by setting Guardian™ Connect to automatically sync to CareLink™ Connect. Our secure online diabetes management platform enables your diabetes care team or your friends and family to have streamlined access to your latest glucose readings, whilst automatic text message alerts can be set up for 5 people giving them added peace of mind, 24/7.

ALWAYS BY YOUR SIDE

With our comprehensive StartRight programme you'll receive all the help you need to as soon as you start. You'll also benefit from our 24/7 MiniMed™ Care Product Support team who are there to answer your questions day or night.



Visit our website for more information and click 'Order Now' for details on how to purchase

guardianconnect.medtronic-diabetes.co.uk or guardianconnect.medtronic-diabetes.ie

Medtronic

United Kingdom

Medtronic UK Ltd.
Building 9, Croxley Park
Watford,
Hertfordshire
WD18 8WW
Diabetes Helpline: +44 (0)1923 212213

Republic of Ireland

Medtronic Ireland Limited
Block 3090-3094
Lake Drive, Citywest Business Campus
Dublin
D24 XN47
Diabetes Helpline: +353 1 511 1444

medtronic-diabetes.co.uk

medtronic-diabetes.ie

Important Safety Information

The Guardian™ Continuous glucose monitoring systems are intended for monitoring interstitial glucose levels in persons with diabetes. These systems provide information regarding potential high and low glucose levels and can provide an indication of when patients should obtain a glucose reading using a home glucose reader. The information provided by continuous glucose monitoring systems is intended to supplement, not replace, readings from a blood glucose meter. Insertion of a glucose sensor may cause bleeding or irritation at the insertion site. Consult a healthcare professional immediately if you experience significant pain or if you suspect that the site is infected. The Guardian™ Connect app requires that you have a functioning mobile electronic device with the correct settings (such as a supported operating system and Bluetooth® enabled) for accurate operation. See website for a list of compatible operating systems.

The patient's device (running the Guardian™ Connect app) must be connected to the Internet via WiFi or a mobile data network (e.g. a 3G or 4G network) to enable a care partner to receive SMS alerts. If the patient receives an alert when the Guardian™ Connect app can't connect to the Internet, the care partner will not be able to receive an SMS alert.

For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use.

Guardian™, CareLink™ are trademarks of Medtronic Inc and its subsidiaries. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google Play and the Google Play logo are trademarks of Google LLC.

UC201710644a - June 2018 ©2018 Medtronic. All rights reserved. Printed in Europe.