

MENU MAP



Home Screen

Press

Suspend Delivery
Audio Options
History
Reservoir & Tubing
Insulin Settings
Sensor Settings
Event Markers
Reminders
Utilities



NAVIGATION

Press from any screen to open the **Menu**.

Press and to scroll through the menu items.

Press on the desired menu item to open.

The scroll bar appears on menus to indicate when additional text is available.

Press to scroll down to view additional items.

Press to scroll back up.

Press to go to previous screen. Hold to return to the **Home** screen.

*Only available when Sensor feature is on.